

Instructions for Use Pre-Trigger Solution

[Product Name]

Pre-Trigger Solution

[Packing Size]

4×125 mL/Pkg, 10×15 mL/Pkg

[Intended Use]

This product is intended to be used in combination with the matching instrument to provide the environment for chemiluminescence reaction.

[Principle]

The pre-trigger solution is added to the reaction mixture to provide an acid environment for releasing the compound marked with acridinium and providing necessary reaction components.

[Main Components]

This product is mainly composed of acid solution containing 1.32% (w/v) hydrogen peroxide.

[Storage Condition and Shelf Life]

Storage condition: Sealed storage at 2~8 °C.

Shelf life: 14 months

After the product is uncapped and placed on the instrument for use, the shelf life is 28 days.

Manufacturing and expiration dates can be found on labels and package.

[Applicable Instrument]

SMART 6500 Automated Chemiluminescence Immunoassay Analyzer

SMART 500S Automated Chemiluminescence Immunoassay Analyzer

Immu F6 Automated Chemiluminescence Immunoassay Analyzer

Immu F6S Automated Chemiluminescence Immunoassay Analyzer

[Product Characteristics]

1. Appearance

(1) The package of the product should be intact without any damage. The appearance is clean and tidy, the label is clear, and no liquid leakage can be found.

(2) The pre-trigger solution is clear liquid without any sediments or floccules.

2. Fill Volume

The fill volume must not be smaller than the value printed on the label.

3. Background RLU

Background RLU of the pre-trigger solution ≤ 500 .

4. Effectiveness

Effectiveness RLU of the pre-trigger solution $\geq 3,000,000$.

[Precautions]

- 1 This product is for IVD use only.
- 2 This product is intended to be used by professional only.
- 3 Do not use the product beyond its expiration date.
- 4 Do not pool the pre-trigger solution from different bottles or different lots.
- 5 Do not violently shake the reagent bottle to avoid bubble generation.
- 6 During use of this product, perform tests by strictly following operation procedures as written in the package insert and guidelines set up by each laboratory.
- 7 Lab operator must wear suitable gloves when using the product. In case of accidental exposure to the reagent, flush the body part with large volume of water immediately. If reagent splashes into eyes, flush eyes with copious of water and consult a doctor immediately.
- 8 All samples and reaction waste must be considered potentially biohazards and be handled in accordance with the local laws and regulations.

[Symbol Description]

	Temperature limit		Date of manufacture
	In vitro diagnostic medical device		Catalogue number
	Batch code		Consult instructions for use
	Use-by date		Authorized representative in the European Community
	This way up		CE marking
	Manufacturer		Unique device identifier
	Importer		

[Basic Information]



MEDCAPTAIN MEDICAL TECHNOLOGY CO., LTD.

12th Floor, Baiwang Research Building, No. 5158 Shahe West Road, Xili, Nanshan, 518055 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA

Telephone: +86-755-26953369

Website: <http://www.medcaptain.com>



Post-sales service: MEDCAPTAIN MEDICAL TECHNOLOGY CO., LTD.

Telephone: +86-755-26953369

Postal code: 518055

Manufacturing Location: Building C, Jiale Science and Technology Industrial Park, Matian Street, Guangming, 518106 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA



R Sight B.V.

Roald Dahllaan 47, 5629 MC, Eindhoven. The Netherlands

[Issued on]

Feb., 25, 2025