KULLANIM AMACI

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Across Gel Forward & Reverse ABO With D^{VI-} / D^{VI+}

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CE
2797

(A / B / AB / D^{VI}- / D^{VI}+ / Ctl / N_{A1} / N_B) | For 'in vitro' diagnostic use only. Read carefully before using the product.

INTENDED PURPOSE

The ABO system was the first human blood group system discovered by Landsteiner⁽¹⁾ and is still the most important system in transfusion practice. Therefore, while classification of ABO system, both forward and reverse grouping must be done. In ABO forward grouping, red blood cells are tested with anti-A, anti-B, anti-AB antibodies, whereas in reverse grouping, ABO antibodies are determined in the serum/plasma using A1 and B red blood cells. Anti D reagent is used for the determination of presence or absence of D (RH1) antigen on the surface of red blood cell.

Blood grouping is life sustaining before blood transfusion. The blood group compatibility between the patient and the donor before the blood transfusion is one of the things that should be done before the blood transfusion. This card is also used for this purpose. It is not used for diagnostic purposes and is not specific to the diagnosis of a specific disease. The ABO system is defined by the presence or absence of the A and/or B antigens on the surface of red blood cells (forward grouping) and by the presence of antibodies in the serum corresponding to the antigen or antigens missing on the surface of red blood cells (reverse grouping).

- This card is suitable for both manual and automatic systems.

- This card group detects parameters qualitatively and semi-quantitatively.

- It is used by professional health teams (medical laboratory technicians, biologists, and other healthcare professionals) working in health institutions such as hospitals and blood centers where blood transfusion treatment is applied.

INDICATION

Across Gel Forward & Reverse ABO with D^{VI}/D^{VI}+ card is used for determination of A (ABO1), B (ABO2), D (RH1) antigens and of possible anti-A and/or anti-B antibodies.

PRINCIPLE

The principle of this test depends on using gel centrifugation technique to detect agglutination reactions which occur when antigens in red blood cells are bound (coated) with corresponding antibodies, present in the reagent, serum and plasma samples.

Across Gel card contains polymerized dextran microspheres (gel) which function as a filter within the buffer solution in all microtubes on the card. Dextran and specific antibody-containing reagent(s) or buffer solution are mixed.

During the centrifugation of the cards, red blood cells agglutinate collected on the gel surface and throughout the length of the column according to their sizes. Non-agglutinated red blood cells go to the bottom of the microtube. The reaction score is determined according to the location of the collapsed cells.

CONTENT

Microtubes of Across Gel Forward&Reverse ABO with D^{VI}/D^{VI}+ card are identified by the front labels of the card:

- **A (ABO1) Microtube:** Monoclonal Anti-A (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15 clone)
- **B (ABO2) Microtube:** Monoclonal Anti-B (mixture of IgM antibodies of murine origin, LB-2, 9621A8 clone)
- **AB (ABO3) Microtube:** Monoclonal Anti-AB (mixture of IgM antibodies of murine origin, BIRMA-1, ES-15, LB-2, 9621A8 clone)
- **D^{VI}(RH1) Microtube:** Monoclonal Anti-D (mixture of IgM antibodies of human origin, RUM-1, P3x61)
- **D^{VI}(RH1) Microtube:** Monoclonal Anti-D (mixture of IgM and IgG antibodies of human origin, clones RUM-1, P3x61, ESD1M, MS-26)
- **Ctl. Microtube:** Neutral gel (buffered solution without antibodies)
- **N Microtubes:** Neutral gel (buffered solution without antibodies)
- All microtubes contain Sodium Azide (NaN₃) as a preservative.
- Reagents are ready for use. Use microtubes immediately when foils are opened.

STORAGE CONDITIONS

Store cards upright at 2-25°C. Do not place cards in front of and/or next to any source of heat (places under direct sunlight, radiators, air conditioners, big devices which radiate heat when operating) or air ducts. Do not freeze cards. When the box is kept in intact package, it is stable until the expiry date specified on the outer label.

OTHER DEVICES INTENDED TO BE USED IN CONJUNCTION WITH THE DEVICE

- Across System semi-automatic or automatic instruments
- Across System Card Centrifuge
- Across System Card Reader
- Across System Dispenser
- Across System Pipette
- Across System Workstation
- Across LISS
- Disposable pipette tips
- Tube for preparation of suspension
- For reverse group test, Across RBC A/B red blood cell reagent

COLLECTION AND PREPARATION OF SAMPLES

Samples

Blood samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used. Do not use haemolysed, cloudy, contaminated blood samples or with clot presence.

Determination of the antigens of ABO/Rh System:

- Fresh blood samples collected in tubes containing anticoagulants like EDTA or citrate can be used.
- Blood samples stored at 2-8°C up to 48 hours after extraction can be used.
- Blood samples stored at 2-8°C, collected from bags containing ACD, CPD, CPDA and SAG-Mannitol which are not beyond the expiry date specified on the label must be used.

Determination of the reverse ABO group:

- Serum or plasma can be used for the determination of Reverse ABO group.
- Serum or plasma samples with anticoagulant (containing plasma) or without anticoagulant (containing serum) can be used.
- Serum or plasma samples stored at 2-8 °C and taken 48 hours ago can be used.
- Serum or plasma samples stored between -20 / -80 °C can be used⁽³⁾.

Sample Preparation

a) Preparing erythrocyte (red blood cell) suspension (To detect ABO/Rh system antigens)

Prepare 1% erythrocyte suspension (Add 5 µl erythrocyte to 0.5 ml Across LISS).

b) Plasma or Serum (Reverse ABO grouping)

Perform necessary centrifuging procedure for clean plasma or serum.

QUALITY CONTROL

It is recommended to include known positive and negative controls in each series of tests. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

TEST PROCEDURE

Across Gel Forward & Reverse ABO with D^{VI}/D^{VI}+ card is used with manual system and Across Auto System semi-automatic or automatic instruments.

Allow the samples and reagents to reach room temperature (18-25°C). Check card condition before use (Do not use the card if you detect bubbles in microtube, absence of supernatant, reduction in gel amount, crack in gel or change in color of the microtube). Across Gel Forward & Reverse ABO with D^{VI}/D^{VI}+ test method is the same principle for manual system and Across Auto System semi-automatic or automatic instruments. For Across Auto System semi-automatic or automatic instruments see the instrument manual.

MANUAL METHOD

- Each card is individually identified by a barcode. Unless a barcode reader is used identify them manually.
- Carefully peel off the aluminium foil that covers the microtubes and keep the card upright to avoid reagent cross contamination.
- Reverse ABO Grouping (N_{A1} - N_B)**
- Homogenise Across Red Blood Cell reagents (A/B) into the N_{A1} / N_B microtubes.
- Pipette 50 µl serum or plasma to microtube N_{A1} / N_B.
- Forward ABO Grouping + D^{VI}/D^{VI}+ Antigen Determination (A-B-AB-D^{VI}-D^{VI}+ -Ctl.)**
- Ensure the resuspension of the red blood cells before use.
- Pipette 50 µl 1% red blood cell suspension into the corresponding microtubes (Avoid contact of pipette tip with the walls or content of microtubes).
- Centrifuge card or cards for 9 mins at 990 rpm.
- Read the results.

RESULTS

Results Reading

Positive	4+	Agglutinated cells in the upper part of the column.
	3+	Medium-sized visible agglutinations in the upper half of the column.
	2+	Small or medium-sized visible agglutinations throughout the column.
	1+	Some small-sized visible agglutinations in the column.
Negative	+/-	Scarce small-sized visible agglutinations in the lower half of the column.
	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

It is recommended that the results are read immediately after centrifuging the cards. Do not leave processed cards in horizontal position. Where necessary, you may store cards and make readings up until 24 hours if cards are kept in vertical position, refrigerated at 2-8 °C and sealed with parafilm or similar material to avoid evaporation of the supernatant.

Interpretation of the results

ABO System

Forward ABO Group				Reverse ABO Group		
A Microtube	B Microtube	AB Microtube	Ctl. Microtube	Microtube N _{A1} + A1 Red Blood Cells Reagent	Microtube N _B + B Red Blood Cells Reagent	ABO Group
0	0	0	0	+	+	0
+	0	+	0	0	+	A
0	+	+	0	+	0	B
+	+	+	0	0	0	AB

Rh System (D antigen)

D ^{VI} - Microtube	D ^{VI} + Microtube	Ctl. Microtube	Interpretation
+	+	0	D positive
0	0	0	D negative
0	+	0	Weak or partial D
+	0	0	

+ = positive

0 = negative

Notes:

- Results are not alone eligible for diagnostic use. They should be evaluated with clinical findings of the patient and other data.
- The meaning of Ctl abbreviation is Control.
- Ctl microtube must be negative. If it is positive, the test is invalid. Wash red blood cells with physiological saline solution, and repeat the process after preparing a new suspension with washed red blood cells. If the repeated control is negative, test results are eligible for evaluation. If it is positive, the test is invalid and must be repeated with a new sample.
- In case of any discrepancy between forward-reverse groups, investigation must be made before giving the samples' test results (patient's age, clinical information, history, etc.).
- In ABO/Rh system : In case of weak positive reactions less than +3, the presence of weak antigens should be investigated.
- D^{VI}/D^{VI}+ microtubes does not react with all D variants so negative reactions with the D^{VI}/D^{VI}+ microtubes must be verified using other reagents and techniques which may detect different variants of the D antigen.
- Any discrepancy between results obtained in D^{VI}- and D^{VI}+ should be interpreted as weak or partial D antigen. Analyzing antigen expression is recommended.

- The O group cells should obtain negative results. Positive results may be caused by unexpected irregular antibodies reacting at room temperature with the corresponding antigens. If positive results are obtained, invalidate the test and proceed to investigate the discrepancy.
- Some individuals with A2 or A2B subgroups may present antibodies anti-A1, where a forward-reverse discrepancy can be seen.

PERFORMANCE CHARACTERISTICS

Specificity and sensitivity of antibodies present in Across Gel Forward & Reverse ABO with D^{VI}/D^{VI}+ card have been studied in a representative number of positive and negative samples.

Antibody	Number of Sample	Sensitivity(I)	Specificity(II)
anti-A	1061	100%	100%
anti-B	1061	100%	100%
anti-AB	1061	100%	100%
anti-D ^{VI} -	1061	98.5%	100%
anti-D ^{VI} +	1061	98.5%	100%

(I) Sensitivity: [(number of true positive results) / (number of true positive results + number of false negative results)] x 100

(II) Specificity: [(number of true negative results) / (number of true negative results + number of false positive results)] x 100

Reverse ABO group

There was compatibility between the results of forward and reverse grouping tests of Across Gel Forward & Reverse ABO with D^{VI}/D^{VI}+ card.

WARNINGS AND PRECAUTIONS

- Using a solution other than Across LISS in preparing red blood cell suspensions may change the reaction.
- Preparing the red blood cell suspension in a concentration other than the specified concentration may result in a false negative test result or a false positive test result.
- Adding volumes other than those specified in the method can alter the reaction.
- Do not use a card with past expiry date.
- Do not use a card whose foil is damaged.
- If dispersed droplets are observed on the microtube due to inappropriate handling or storage, it is recommended that cards are centrifuged before use. If droplets do not precipitate, do not use the card.
- Abnormal concentrations of serum proteins, presence of macromolecular solutions in serum or plasma, or presence of Wharton jelly in umbilical cord blood samples may result in non-specific agglutination of red blood cells. In such cases, it is recommended that red blood cells are washed before the test⁽²⁾.
- Double population may be seen in patients who have undergone blood transfusion or bone marrow transplantation⁽²⁾.
- As red blood cells may be completely coated with antibodies in patients with high cold antibody effect, spontaneous agglutination may occur⁽²⁾.
- All blood products and samples should be handled carefully as if they have potential for transmitting of disease.
- The product is recommended to be disposed of into a special container used for biological waste after use.
- If you encounter any adverse event or need further information regarding the use of this product, consult the authorized distributor in your country.
- Report any serious adverse events that occur with the device to the manufacturer (diapro@diapro.com.tr / +90(262)751-53-34) and the competent authority in which you reside.

REFERENCES

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 10th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005.
3. WHO, guidelines and principles for safe blood transfusion practice module 3.

PACKAGING

Across Gel Forward & Reverse ABO With D^{VI}- / D^{VI}+

(A / B / AB / D^{VI}- / D^{VI}+ / Ctl / N_{A1} / N_B)

Ref No: 810201 ---- 50 Cards

Document No: TK3539

This document is available in several languages. Translations have been made from the IFU in Turkish. In case of doubt or discrepancy, statement in master Turkish document shall be taken into consideration.

<https://ec.europa.eu/tools/eudamed>

SSP Ref No: SSP-AGC-810201

BUDI:868019639TDAG1DNU

IVD

"In Vitro" Diagnostic Medical Device

LOT

Batch Code

Use by date

Temperature Limit

Manufacturer

Consult Instructions for Use

REF

Catalogue Number

Do not re-use

Do not use if package is damaged and consult

Contains biological material of animal origin

Contains biological material of human origin