




HS-OT10J

Specification

- Direct illumination
- 3X magnification for the viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- The clip with integrated on/off switch
- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- Colors to choose
- 直接照明
- 3倍放大观察窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 固定夹扣式on/off 开关
- 3200k和6000k色温LED灯和XHL 卤素灯可选
- 五种颜色可选择



Direct illumination

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
193x102x32	400x215x195	20	4.9/4.5

Gross and net weight(Kgs/Ctn)Battery containing

Specification

- Direct illumination
- 3X magnification for the viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- Standard handle,handy and convenient
- The clip with integrated on/off switch
- 直接照明
- 3倍放大观察窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 3200k和6000k色温LED灯和XHL 卤素灯可选
- 标准型手柄,轻巧方便
- 固定夹扣式on/off 开关



Direct illumination



HS-OT10G

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
190x110x40	400x240x230	20	6.8/6.3



Gross and net weight(Kgs/Ctn)Battery not included



HS-OT10N

Specification

- Fiber optic illumination
- 3X magnification for the viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- Supplied by 2xAA alkaline batteries
- 光纤照明
- 3倍放大观察窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 3200k和6000k色温LED灯和XHL 卤素灯可选
- 2×AA碱性电池



Fiber optic illumination


Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
193x102x32	400x215x195	20	4.9/4.5

Gross and net weight(Kgs/Ctn)Battery containing

Specification

- Fiber optic illumination
- 3X magnification for the large size viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- The socket of ear insufflation bulb
- 光纤照明
- 3倍放大的大观察窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 3200k与6000k色温LED灯和XHL 卤素灯可选
- 气囊球插口



Fiber optic illumination



HS-OT10S

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
193x102x32	400x215x195	20	4.5/4.1

Gross and net weight(Kgs/Ctn)Battery not included