

NEW + LED

HS-OTP10L-EY



- A deluxe LCD screen handle
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- A head of fiber optic otoscope
- A head of direct ophthalmoscope
- A head of nasal speculum
- A head of tongue depressor holder
- 5 pcs tongue depressor

HS-OTP10L-KE



- A deluxe LCD screen handle
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- A head of fiber optic otoscope
- A head of direct ophthalmoscope
- A head of nasal speculum
- A head of laryngeal mirror and reflector lens

HS-OTP10L-KY



- A deluxe LCD screen handle
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- A head of fiber optic otoscope
- A head of direct ophthalmoscope
- A head of tongue depressor holder
- A head of laryngeal mirror and reflector lens
- 5 pcs tongue depressor

HS-OT10L-KEY



- A deluxe LCD screen handle
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- A head of fiber optic otoscope
- A head of tongue depressor holder
- A head of laryngeal mirror and reflector lens
- 5 pcs tongue depressor

- The handle of the above series of products can be selected as a standard type of handle
- 上面整套系列产品的柄可选配为标准手柄



Package content

Model	Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
HS-OTP10L-EY	275x160x56	340x255x310	10	9.2/8.5
HS-OTP10L-KE	275x160x56	340x255x310	10	8.8/8.1
HS-OTP10L-KY	275x160x56	340x255x310	10	8.2/7.5
HS-OT10L-KEY	275x160x56	340x255x310	10	9.5/8.7

Gross and net weight: Kgs(Ctn)Battery not included



HS-OT10C

- A rechargeable fiber optic otoscope  
一个可充电光纤检耳镜
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups  
一组2.2mm/3mm/4mm/5mm/9mm耳套
- AA Ni-Mh2.4V rechargeable battery  
2.4V AA镍氢充电电池
- Single-charger  
单孔充电器

- Search bar and lens are removable  
可拆卸的探索杆和镜头
- Straight rod or Curved rod are optional  
弯杆直杆可选择



HS-OP10C

- 5 different apertures  
五种光阑
- 19 kinds of adjustable corrective lenses,  
Range of lenses: -20d ~ +20d  
19种屈光镜片, 范围: -20D ~ +20D

- 2x optical glass observation window lens  
2x光学玻璃观察窗镜片
- It can be upgraded to the HS-OT10H  
可升级为HS-OT10H手术耳镜
- 5 pcs 9mm nose tips  
5个9mm窥鼻器



HS-OT10K3-C

- 5 different apertures  
五种光阑
- 19 kinds of adjustable corrective lenses,  
Range of lenses: -20d ~ +20d  
19种屈光镜片, 范围: -20D ~ +20D



HS-OT10E3-C

- Patented tongue-depressor ejectable device  
拥有专利的压舌板弹出装置

Package content

Model	Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
HS-OT10K3-C	206x173x88	460x300x330	12	11.6/10.6
HS-OT10C	206x173x88	460x300x330	12	11.5/10.5
HS-OP10C	206x173x88	460x300x330	12	11.0/10.0
HS-OT10E3-C	206x173x88	460x300x330	12	11.8/10.8
HS-OT10E3-C	206x173x88	460x300x330	12	11.9/10.9

Gross and net weight: (Kgs/Ctn)Battery containing

- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to options  
3200k和6000k色温LED灯和XHL卤素灯可选