

BESDATA® M7 7inch LCD Monitor



HDMI output

16GB Storage Card

Compatible with BESDATA Reusable /
Disposable laryngoscopes, Video Stylet,
Double lumen tube etc.

Technical Data Sheet

Parts	Description	Technical Index
Monitor	Display	7 inch LCD
	Resolution Ratio	1280*720(RGB)
	Viewing angle	>80°
	Power	<2W
	Weight	750g
Camera	Illuminance	>150 LUX
	Resolution Ratio	1280*720 px
Image/Video function	Image/ Video function	Yes
	Output	SD Memory Card And HDMI Video Output
		Easy to storage and establish image/video
	Memory type	Micro SD card 16GB
	Video/Photo Play back	Yes
Batterty	Battery Type	Rechargeable lithium battery
	Capacity	4300mAh
	Cycle life of battery	>300 times
	Battery working time	>120min
	Charging time	<150min