INT-BRO-202411-20

Anbio PRODUCT BROCHURE

ANBIO BIOTECHNOLOGY

The healthcare landscape is evolving.
We advance diagnostics by providing accessible diagnostic solutions at the forefront of science.

DIAGNOSTICS BEYOND

FIA ChLIA LFIA ABO&RhD Dry Chemistry Molecular Dx Vet Dx

CONTACT US

Contact us at www.anbio.com/contact to learn more or visit www.anbio.com

Subscribe us to get latest updates of Anbio

©2024 Anbio Biotechnology. All rights reserved.


ABOUT ANBIO

Breaking Barriers in Diagnostics

Anbio Biotechnology (Anbio) is a leading global medical device company with a focus on in vitro diagnostics. Our goal is to revolutionize the diagnostics market by personalizing and decentralizing current diagnostic solutions, enabling rapid diagnosis and improving patient prognosis. We achieve this by offering cutting-edge laboratory, wellness, at-home, and point-of-care (POCT) in vitro diagnostic (IVD) solutions that are rapid, accessible, and affordable on a global scale.

Since its establishment, Anbio has been committed to delivering the most advanced diagnostic solutions while ensuring cost-effectiveness for our customers. Our extensive capabilities in research and development, manufacturing, and global product distribution logistics allowing us to continuously focus on cost reduction while offering a wide range of mature laboratory diagnostics and point-of-care technology (POCT) products. Our product range spans across cancer, cardiovascular conditions, pharmacogenomics, hormones, infectious diseases, inflammation, and drugs of abuse tests.

ANBIO'S GLOBAL PRESENCE



AFRICA

- Kenya
- Nigeria
- Egypt
- Uganda
- Ghana
- Ethiopia
- .....

AMERICA

- USA
- British Virgin Islands
- Cayman Islands
- Bolivia
- Brasil
- Colombia
- .....

EUROPE

- United Kingdom
- France
- Germany
- Italy
- Spain
- Portugal
- .....

APAC

- Hong Kong SAR
- Mainland China
- Malaysia
- Indonesia
- Philippines
- Vietnam
- Thailand
- .....

Global IVD Supply and Manufacturing Partner

- Over 300 IVD Products
- High Manufacturing Capability
- Over 1 Million Sqft (100,000 m^2) Manufacturing Site
- Quality Management System
  ISO13485
  MDSAP

- Technology Licensing and Sharing with Global Partners
- Anbio can deliver to Global Partners, monthly:

  700 million Lateral Flow Tests
  700 million FIA Tests
  500 million Chemiluminescent Tests
  300 million Molecular Tests
  100 million Molecular POCT Tests



Immunoassay Diagnostic Platform

- **Multiple options**
  Suitable for multiple applications
- **Cost Effective**
- **Accurate**

Chemiluminescence Immunoassay



Lateral Flow Immunoassay



Fluorescent Immunochromatography Assay



Anbio’s FIA Solution

The Anbio FIA Platform is a handheld, easy-to-use, affordable, and rapid point-of-care FIA that utilizes an LED light source from the reader for excitation of fluorescent dye for immunochromatographic qualitative or quantitative testing of analytes in whole blood and urine samples, including hormones detection, myocardial disease detection, infectious disease detection, and tumor-related antigens detection.

FLUORESCENCE IMMUNOASSAY DIAGNOSTIC SERIES



AF-100 Single Channel Benchtop FIA Analyzer



- Single Channel
- Automatically Identify Test Item and Disposal
- Friendly Touch Screen and Intuitive Interface
- Easy to Operate
- Deliver Rapid and Accurate Results in 3-15 Mins

Model	AF-100
Test speed	240 tests/hour
Discard card mode	Automatic Discard
Test items	Automatic identification, supporting more than 80 items
Memory capacity	8000 results (customizable)
Display screen	7-inch color screen
Control mode	Capacitor Touch
Printing	Built-in thermal printer
Data transmission	USB and RS-232
Size	265*240*130mm

AF-100S Handheld FIA Analyzer



- Single Channel
- Palm Size with Built-In Battery, Enable Medical Testing Anywhere, Anytime
- Friendly Touch Screen and Intuitive Interface
- Easy to Operate
- Compatible with Bluetooth, Transfer Data Wireless

Model	AF-100S
Test speed	240 tests/hour
Test items	Automatic identification, supporting more than 80 items detection
aMemory capacity	8000 results (customizable)
Display screen	4.3-inch touchable True color display 800*480 pixel
Control mode	Capacitor Touch
Printing	Adapted Bluetooth Thermal Printer
Data transmission	USB Type C
Built-in battery	7.4V/2800mAh Continuous uninterrupted testing>800tests Standby 20 days
Charge time	<120mins
Size	195*100*70mm

AF-1200 Multi-channel Benchtop FIA Analyzer



- 12 Testing Channels, Enhance Testing Efficiency
- Easy to Operate
- Friendly Touch Screen and Intuitive Interface
- Deliver Accurate and Reliable Result in 3 -15 Mins

Model	AF-1200
Test speed	48-240 tests/hour
Incubation Temperature	26℃
Test Method	Timed incubation for automatic detection, simultaneously supporting 12 different test items.
Test items	Automatic Identification, supporting more than 80 items
Memory capacity	8000 results (customizable)
Display screen	10.1-inch color screen
Control mode	Capacitor Touch
Printing	Built-in thermal printer
Data transmission	USB and RS-232
Size	425*260*175mm

AF-100C Single-Channel Benchtop FIA Analyzer



- Single Channel
- Built-In Battery, Enable Medical Testing Anywhere, Anytime
- Friendly Touch Screen and Intuitive Interface
- Easy to Operate

Model	AF-100C
Test speed	180 tests/hour
Test items	Multiple specimen Type Detection, Automatic Identification, supporting more than 80 items
Memory capacity	8000 results
Display screen	7' IPS Capacilive touch screen
Control mode	Capacitor Touch
Printing	Built-in thermal printer
Data transmission	Onboard screen / Printer / USB
Results	Results in 3-15 Minutes
Net Weight	1.2Kg
Size	250*120*90mm

Fluorescence Immunoassay kit series



Anbio Fluorescence Immunoassay Series - Ensuring Health Anywhere Anytime

- Unique Design for Portable Diagnostic Testing
- Rapid and Accurate Results in 15 Minutes
- Automatic Cartridge Identification and Disposal
- Compatibility with Serum, Plasma, Whole Blood, Urine, Feces and Other Sample Extracts
- Ideal for Mobile Testing in Laboratory-Scarce Areas
- Suitable for Near-Patients Testing, Ambulances, and Mobile Healthcare Platforms
- User-Friendly Touchscreen Graphic User Interface


FIA Parameters

REF NO	Marker	Product Name	Type	Range of Linearity	Specimen Type	Specimen Volume	Reaction Time	Qualification
A02084	Malaria (Pf/Pan)	Malaria (Pf/Pan) Rapid Test Kit	Parasitic Infections	Qualitative detection	WB	5μL	15mins	CE
A02083	Malaria (Pf/Pv)	Malaria (Pf/Pv) Rapid Test Kit		Qualitative detection	WB	5μL	15mins	CE
A02085	Dengue IgM/IgG	Dengue IgM/IgG Rapid Test Kit	Viral Infections	Qualitative detection	S/P/WB	10μL/10μL/10μL	15mins	CE
A02107	Dengue NS1	Dengue NS1 Rapid Test Kit		Qualitative detection	S/P/WB	50μL/50μL/50μL	15mins	CE
A02034	Hepatitis B surface antigen	HBsAg Rapid Test Kit		Qualitative detection	S/P/WB	50μL/50μL/50μL	15mins	
A02169	Hepatitis B surface antibody	HBsAb Rapid Test Kit		5-1000mIU/mL	S/P/WB	50μL/50μL/50μL	15mins	
A02035	Hepatitis C Virus	HCV Rapid Test Kit		Qualitative detection	S/P/WB	50μL/50μL/50μL	15mins	
A06021501	COVID-Neutralizing Antibody	Neutralizing Antibody of SARS-CoV-2 Test(Fluorescence Immunochromatographic Assay)		10.00 ~1000.00 IU/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A6061293	COVID-Ag	Rapid COVID-19 Antigen Test		Qualitative detection	NP/OP	Refer to IFU	15mins	
A02105	Influenza A/B	Influenza A/B Rapid Test Kit		Qualitative detection	NP/OP	Refer to IFU	15mins	CE
A02108	Respiratory Syncytial Virus Ag	RSV Ag Rapid Test Kit		Qualitative detection	NP/OP	Refer to IFU	15mins	CE
A02099	Human immunodeficiency virus Ab	Anti-HIV (1&2) Rapid Test Kit		Qualitative detection	S/P/WB	50μL/50μL/50μL	15mins	
A02165	Human immunodeficiency virus Ag/Ab	HIV Ag/Ab Rapid Test Kit		Qualitative detection	S/P/WB	50μL/50μL/50μL	15mins	
A02115	Treponema Pallidum Ab (Syphilis)	Anti-TP Rapid Test Kit		Qualitative detection	S/P/WB	50μL/50μL/50μL	15mins	CE
A02082	Helicobacter pylori Ab	Anti-HP Rapid Test Kit	Bacterial Infections	Qualitative detection	S/P/WB	50μL/50μL/50μL	15mins	CE
A02081	Helicobacter pylori Ag	HP Rapid Test Kit		Qualitative detection	Stool	100ul/50mg	15mins	CE
A02167	Group A Streptococcus	Group A Streptococcus Test Kit		Qualitative detection	Throat Swab	Refer to IFU	5mins	CE
A02109	Fecal Occult Blood	FOB Rapid Test Kit	Degistive Tract Disease	0.05 ~ 20.00 μg/mL	Stool	100ul/50mg	15mins	CE
A02086	Calprotectin	Calprotectin Rapid Test Kit		10.00 ~ 1000.00 μg/g	Stool	100ul/50mg	15mins	CE
A02103	Troponin I Ultra	Troponin I Ultra Rapid Test Kit	Cardiac Marker	0.03 ~ 20.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02001	N-terminal pro-B-type natriuretic peptide	NT-proBNP Rapid Test Kit		150.00 ~ 35000.00pg/mL	S/P/WB	50μL/50μL/50μL	15mins	CE/NMPA
A02003	Troponin I	cTnI Rapid Test Kit		0.10 ~ 50.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE/NMPA
A02004	Troponin T	cTnT Rapid Test Kit		50.00 ~ 10000.00 pg/mL	S/P/WB	50μL/50μL/50μL	15mins	CE/NMPA
A02005	Myoglobin	Myo Rapid Test Kit		10.00 ~ 400.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02006	Creatine Kinase-MB	CK-MB Rapid Test Kit		2.00 ~ 100.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02010	Soluble suppression of tumorigenicity 2	sST2 Rapid Test Kit		5.00 ~ 400.00ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE/NMPA
A02046	Heart-type Fatty Acid Binding Protein	H-FABP Rapid Test Kit		1.00 ~ 120.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02002	CK-MB/cTnI/Myo	CK-MB/cTnI/Myo Rapid Test Kit		CK-MB: 2.00 ~100.00 ng/mL; cTnI: 0.10 ~ 50.00 ng/mL; Myo: 10.00 ~ 400.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE/NMPA
A02008	D-Dimer	D-Dimer Rapid Test Kit		0.10 ~ 10.00 mg/L FEU	P/WB	50μL/50μL	15mins	CE/NMPA
A02016	S100 calcium-binding protein beta	S100-β Rapid Test Kit	Cerebral Apoplexy	0.08 ~ 10.00ng/mL	S	50μL	15mins	CE
A02011	C-reactive protein	CRP Rapid Test Kit	Inflammation	0.50 ~ 200.00 mg/L	S/P/WB	10μL/10μL/10μL	3mins	CE/NMPA
A02012	Procalcitonin	PCT Rapid Test Kit		0.10 ~ 100.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE/NMPA
A02013	Serum amyloid A	SAA Rapid Test Kit		3.00 ~ 200.00 mg/L	S/P/WB	10μL/10μL/10μL	3mins	CE/NMPA
A02014	Interleukin-6	IL-6 Rapid Test Kit		3.00 ~ 4000.00 pg/mL	S/P/WB	50μL/50μL/50μL	15mins	CE/NMPA
A02015	Heparin-Binding Protein	HBP Rapid Test Kit		5.00~300.00ng/mL	P(only EDTA)	50μL	15mins	CE/NMPA
A02064	Procalcitonin/Serum amyloid A	CRP/SAA Rapid Test Kit		CRP: 0.50 ~ 200.00 mg/L SAA: 3.00 ~ 200.00 mg/L	S/P/WB	10μL/10μL/10μL	3mins	CE
A02049	Free Triiodothyronine	FT3 Rapid Test Kit	Thyroid Hormones	1.00 ~ 50.00 pmol/L	S/P	50μL/50μL	15mins	CE
A02050	Free thyroxine	FT4 Rapid Test Kit		1.00 ~100.00 pmol/L	S/P	50μL/50μL	15mins	CE
A02022	Triiodothyronine	TT3 Rapid Test Kit		0.80 ~ 10.00 nmol/L	S/P	50μL/50μL	15mins	CE
A02023	Thyroxine	TT4 Rapid Test Kit		15.00 ~ 320.00 nmol/L	S/P	50μL/50μL	15mins	CE
A02024	Thyroid-stimulating hormone	TSH Rapid Test Kit		0.10 ~100.00 mIU/L	S/P/WB	50μL/50μL/50μL	15mins	CE
A02027	Follicle-Stimulating Hormone	FSH Rapid Test Kit	Reproductive Hormones	1.00 ~ 200.00 mIU/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02042	Progesterone	PROG Rapid Test Kit		1.00 ~ 60.00 ng/mL	S/P	50μL/50μL	15mins	CE
A02043	Testosterone	TES Rapid Test Kit		0.20 ~16.00 ng/mL	S/P	50μL/50μL	15mins	CE

FIA Parameters

REF NO	Marker	Product Name	Type	Range of Linearity	Specimen Type	Specimen Volume	Reaction Time	Qualification
A02026	Luteinizing Hormone	LH Rapid Test Kit	Reproductive Hormones	1.00 ~ 200.00 mIU/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02028	Prolactin	PRL Rapid Test Kit		1.00 ~ 200.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02007	Estradiol	E2 Rapid Test Kit		30.00 ~ 3000.00 pg/mL	S/P	50μL/50μL	15mins	CE
A02025	β-human chorionic gonadotrophin	β-HCG Rapid Test Kit		2.00 ~ 200,000.00 mIU/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02041	Anti-Müllerian hormone	AMH Rapid Test Kit		0.20 ~ 20.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02018	Cortisol	Cor Rapid Test Kit	Adrenal Hormones	5.00 ~ 2000.00 nmol/L	S/P/WB	50μL/50μL/50μL	15mins	CE
A02032	Carcinoembryonic antigen	CEA Rapid Test Kit	Cancer Marker	1.00 ~ 500.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02051	Carbohydrate antigen 125	CA125 Rapid Test Kit		2.00 ~ 600.00 U/ml	S/P/WB	50μL/50μL/50μL	15mins	CE
A02052	Carbohydrate antigen 15-3	CA15-3 Rapid Test Kit		3.00 ~ 300.00 U/ml	S/P/WB	10μL/10μL/10μL	15mins	CE
A02053	Carbohydrate antigen 19-9	CA19-9 Rapid Test Kit		2.00 ~ 700.00 U/ml	S/P/WB	50μL/50μL/50μL	15mins	CE
A02048	Alpha-Feto Protein	AFP Rapid Test Kit		5.00 ~ 500.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02130	Pepsinogen I	PGI Rapid Test Kit		2.00 ~ 200.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02131	Pepsinogen II	PGII Rapid Test Kit		1.00 ~ 100.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02020	Gastrin-17	G-17 Rapid Test Kit		1.00 ~ 40.00 pmol/L	S/P/WB	50μL/50μL/50μL	15mins	CE/NMPA
A02021	Pepsinogen I/ II	PGI/PGII Rapid Test Kit		PGI: 2.00 ~200.00 ng/mL; PGII: 1.00 ~100.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE/NMPA
A02030	Prostate-specific antigen	PSA Test Kit		0.10 ~ 100.00 ng/ml	S/P/WB	50μL/50μL/50μL	15mins	
A02031	free Prostate-specific antigen	fPSA Test Kit		0.1 ~ 60.00 ng/ml	S/P/WB	50μL/50μL/50μL	15mins	
A02017	Glycated hemoglobin	HbA1c Rapid Test Kit	Diabetes	3.00% ~ 15.00%	WB	10μL	5mins	CE
A02078	Insulin	INS Rapid Test Kit		0.50 ~ 500.00 μU/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02033	Immunoglobulin E	IgE Rapid Test Kit	Immunological Function	1.00 ~1000.00 IU/mL	S/P/WB	10μL/10μL/10μL	15mins	CE
A02065	Cyclic Citrullinated Peptide antibody	Anti-CCP Rapid Test Kit	Rheumatoid Arthritis	10.00 ~ 500.00 U/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02089	Rheumatoid factor	RF Rapid Test Kit		10.00 ~ 160.00 IU/mL	S/P/WB	10μL/10μL/10μL	5mins	CE
A02166	Anti-Streptolysin O	ASO Rapid Test Kit		20.00 ~ 500.00 IU/mL	S/P/WB	10μL/10μL/10μL	5mins	CE
A02088	Insulin-like Growth Factor 1	IGF-1 Rapid Test Kit	Growth Factor	15.00 ~ 1000.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02019	25-hydroxy Vitamin D	25-(OH) VD Rapid Test Kit	Bone Matabolism	10.00 ~100.00ng/mL	S/P/WB	10μL/10μL/10μL	15mins	CE/NMPA
A02047	Neutrophil gelatinase-associated lipocalin	NGAL Rapid Test Kit	Renal Function	10.00 ~ 1500.00 ng/mL	Urine	50μL	15mins	CE
A02040	Microalbuminuria	MAU Rapid Test Kit		5.00 ~ 200.00 mg/L	Urine	50μL	15mins	CE
A02029	Ferritin	Ferritin Rapid Test Kit	Anemia	2.00 ~1000.00 ng/mL	S/P/WB	50μL/50μL/50μL	15mins	CE
A02079	Vitamin B12	VB12 Rapid Test Kit		100.00 ~ 2000.00 pg/mL	S/P	50μL/50μL	10mins	CE
A02080	Folate	Folate Rapid Test Kit		1.00 ~ 30.0 ng/mL	S	10μL	15mins	CE
A02054	Tetrahydrocannabinol	THC Rapid Test Kit	DOA-hair	0.02 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02055	Caffeine	CAF Rapid Test Kit		1.00 ~ 50.00 ng/mg	Hair	5-8mg	5mins	CE
A02056	Methcathinone	MCAT Rapid Test Kit		0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02057	Tramadol	TRA Rapid Test Kit		0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02058	3,4-methylenedioxy-n-methylamphetamine	MDMA Rapid Test Kit		0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02061	Synthetic Cannabinoids	Synthetic Cannabinoids Rapid Test Kit		0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02106	Cocaine	COC Rapid Test Kit		0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02036	Methamphetamine	MET Rapid Test Kit		0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02037	Morphine	MOP Rapid Test Kit		0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02038	Ketamine	KET Rapid Test Kit		0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02059	Methamphetamine/Ketamine	Drug Abuse Rapid Test Kit		MET: 0.10 ~ 10.00 ng/mg MOP: 0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02039	Methamphetamine/Morphine/ Ketamine	Drug Abuse Rapid Test Kit		MET: 0.10 ~ 10.00 ng/mg MOP: 0.10 ~ 10.00 ng/mg KET: 0.10 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02077	COC/MDMA/THC	Drug Abuse Rapid Test Kit		COC: 0.10 ~ 10.00 ng/mg MDMA: 0.10 ~ 10.00 ng/mg THC: 0.02 ~ 10.00 ng/mg	Hair	5-8mg	5mins	CE
A02168	Etomidate	Etomidate Rapid Test Kit		0.1-10.0 ng/mg	Hair	5-8mg	5mins	CE

Immunofluorescence Quality Control Series



No	Ref No	Product	Methodology	Specification	Form
1	0702002	CK-MB/cTnI/Myo Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
2	0702001	NT-proBNP Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
3	0702003	cTnI Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
4	0702103	Troponin I Ultra Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
5	0702004	cTnT Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
6	0702005	Myo Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
7	0702006	CK-MB Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
8	0702008	D-Dimer Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
9	0702022	TT3 Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
10	0702023	TT4 Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
11	0702024	TSH Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
12	0702049	FT3 Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
13	0702050	FT4 Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
14	0702026	LH Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
15	0702027	FSH Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
16	0702042	PROG Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
17	0702043	TES Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
18	0702025	β-HCG Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
19	0702041	AMH Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
20	0702011	CRP Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
21	0702012	PCT Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
22	0702013	SAA Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
23	0702014	IL-6 Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
24	0702032	CEA Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
25	0702048	AFP Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
26	0702030	PSA Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
27	0702031	fPSA Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
28	0702051	CA125 Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
29	0702052	CA15-3 Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
30	0702053	CA19-9 Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
31	0702065	Anti-CCP Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
32	0702089	RF Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
33	0702166	ASO Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
34	0702033	IgE Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
35	0702017	HbA1c Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
36	0702029	Ferritin Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
37	0702019	25- (OH) VD Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
38	0702079	VB12 Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
39	0702078	INS Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
40	0702040	MAU Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
41	0702109	FOB Quality Control	Quantitative detection	0.6ml*2/box	Lyophilized microspheres
42	0702167	Group A Streptococcus Quality Control	Qualitative detection	0.6ml*2/box	Lyophilized microspheres
43	0702107	Dengue NS1 Quality Control	Qualitative detection	0.6ml*2/box	Lyophilized microspheres
44	0702034	HBsAg Quality Control	Qualitative detection	0.6ml*2/box	Lyophilized microspheres
45	0702035	HCV Quality Control	Qualitative detection	0.6ml*2/box	Lyophilized microspheres
46	0702105	Influenza A/B Quality Control	Qualitative detection	0.6ml*2/box	Lyophilized microspheres
47	0702108	RSV Ag Quality Control	Qualitative detection	0.6ml*2/box	Lyophilized microspheres
48	0702081	HP Ag Quality Control	Qualitative detection	0.6ml*2/box	Lyophilized microspheres