Rainhome Medical



Wounds repair specialist |
Focus on wounds repair for 20 years

About Rainhome

Guangzhou Rainhome Pharm & Tech Co., Ltd (Rainhome Medical) established in 2004 , focusing on wounds repair for 20 years , is a national hi-tech manufacturer specializing in design, development, production and sales of medical wound dressings and disinfectants products.

Our factory is located in Development Economic District of Guangzhou city, having been building in line with the management norms of medical device manufacturing enterprises. It has 100000 level purification production workshop, production line and 10000 level Microbial Laboratory, and we have approved CE, ISO EN13485 and passed the quality system assessment of medical device manufacturers.



Rainhome medical has domestic leading scientific research strength and many independent intellectual property rights, and obtained the intellectual property management system certification in 2017. So far we have approved complete certificates such as FDA, CE, ISO EN13485, MSDS, Medical Device Production License and so on. Our main products include medical function dressing, negative pressure wound therapy device, NPWT dressing kit, Alginate dressing, Chitosan dressing, Foam dressing, Hydrocolloid dressing, Gelling fiber dressing, Hydrocolloid oil gauze, medical scar paster, Saline Nasal Spray, Ultrasound Gel, Hand Sanitizer and so on, we have been exporting the goods to different place of the world.



润虹医药

Certificates

We have got ANVISA, SFDA, MDA, INVIMA Certificates and others registrations in others markets.


























Cooperated with wound care hospitals

NPWT system for wound care healing



Cooperated with wound care center










NPWT system therapy clinical case



Catalog

Rainhome is a hi-tech factory specializing in design, development, production and sales of medical device and dressing, as well as disinfectants and skin care.



Welsuc Negative Pressure Wound Therapy Device



RH-S01 NPWT



RH-S01I NPWT



RH-450 NPWT



RH-450I NPWT



RH-1000 NPWT

- Instrument classification：Class II
- Scope of Application ：suitable for Acute & Chronic wounds, diabetic foot ulcers, venous ulcers, bedsores, neuropathic ulcers, donor sites, omphalocele, post-operative abdominoplasty , etc.

Product name	Treatment model	Model	Power source	Air flow rate	Pressure range	Normal working noise	Shelf life	Dimensions (cm)	Weight(g)	Canister capacity	Compatible Dressings' sizes	Device
NPWT system	Continunous therapy	RH-S10	Internal Battery	<2L/min	-80 or -100 or -125mmHg(-6.7~-29.9kPa)	≤45dB	3 years	13.9×10.4×11.6	92	Without canister but can be customized	10x20cm, 15x15cm,15x20cm,15x30cm, 17x22cm,10x40cm	Disposable
NPWT system	Continunous/Dynamic/ Intermittent therapy	RH-S10I	Recharge battery	<2L/min	-40~-300mmHg(-6.7~-29.9kPa)	≤52dB	3 years	8.5×7.5×2	260	50ML, 80ML	5x5x3cm, 10x7x3cm, 15x10x3cm, 20x10x3cm, 25x15x3cm, 30x20x3cm, others customized	Reusable
NPWT device	Continunous/Dynamic/ Intermittent therapy	RH-450	100~240VAC & Recharge battery	<6L/min	-40~-300mmHg(-6.7~-29.9kPa)	≤45dB	3 years	13.9×10.4×11.6	720	140ml, 500ml, 1000ml	5x5x3cm, 10x7x3cm, 15x10x3cm, 20x10x3cm, 25x15x3cm, 30x20x3cm, others customized	Reusable
NPWT device	Continunous/Dynamic/ Intermittent therapy	RH-1000	100~240VAC & Recharge battery	<20L/min	-40~-300mmHg(-6.7~-66.5kPa)	≤55dB	3 years	25×20×25	1500	1000ml, 2000ml	5x5x3cm, 10x7x3cm, 15x10x3cm, 20x10x3cm, 25x15x3cm, 30x20x3cm, others customized	Reusable
NPWT system with instillation	Five treatment models	RH-450I	100~240VAC & Recharge battery	<10L/min	-40~-300mmHg(-6.7~-29.9kPa)	≤45dB	3 years	13.9×10.4×11.6	920	140ml, 500ml, 1000ml	5x5x3cm, 10x7x3cm, 15x10x3cm, 20x10x3cm, 25x15x3cm, 30x20x3cm, others customized	Reusable

Indications for Use:The welsuc Negative Pressure Wound Therapy (NPWT) System is indicated for patients who would benefit from a suction device (Negative Pressure Wound Therapy) as when used on open wounds it creates an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudates and infectious material.

RH 润虹医药 RAINHOME

Welsuc NPWT Device with WIFI / 4G

It can be managed all of the machines via computor and phone.



RH-450 NPWT




RH-1000 NPWT



RH-450I NPWT

- Instrument classification：Class II
- Scope of Application ：suitable for Acute & Chronic wounds, diabetic foot ulcers, venous ulcers, bedsores, neuropathic ulcers, donor sites, omphalocele, post-operative abdominoplasty , etc.



The WIFI or 4G function will be including management of each machine status, therapy situation, alarm, location of each machine and so on.

Difference: if machines with WIFI system, it no need to insert the phone card or internet card into the machines, but it could not check the location. If with 4G, need phone or internet card, and locationing each machine.

RH RAINHOME 润虹医药

Welsuc Negative Pressure Wound Therapy Device

Clinical cases:

Application for the wounds as below:

- Acute or subacute wounds.
- Chronic wounds
- Dehisced wounds
- Pressure ulcers
- Diabetic/Neuropathic ulcers
- Venous insufficiency ulcers
- Traumatic wounds
- Partial thickness burns
- Flaps and grafts
- Closed surgical incisions















- Instrument classification：Class IIa for the device
- Negative pressure wound therapy (NPWT), also referred to as vacuum-assisted wound closure, is a therapeutic dressing system in which negative pressure is continuously or intermittently applied to the surface of a wound. The system includes dressings, a suction pump, tubing and a collection chamber.


Accessories of RH-450 & RH-450I NPWT device





Canister 140mL & 500mL & 1000ML

Just for RH-450 & RH-450I NPWT device



Drainage tubing

Just for RH-450I device

- Instrument classification：Class I
- Scope of Application ：suitable for Acute & Chronic wounds, diabetic foot ulcers, venous ulcers, bedsores, neuropathic ulcers, donor sites, omphalocele, post-operative abdominoplasty , etc.
- **It is combined with welsuc NPWT Device.**

Product name	Model/Size	Fuctures	Advantages	Shelf life	Dimensions (cm)	Weight(g)	Packaging Specifications
Canister	GF-1/140ml	Disposable, Un-sterile (if need sterile, it's available)	With charcoal and membrane filters and solidifier	3 years	4.1×9.5×10.6	127	1pcs/bag
Canister	GF-1/500ml	Disposable, Un-sterile (if need sterile, it's available)	With charcoal and membrane filters and solidifier	3 years	9.8×9.5×10.6	185	1pcs/bag
Canister	GF-1/1000ml	Disposable, Un-sterile (if need sterile, it's available)	With charcoal and membrane filters and solidifier	3 years	18.5×15.5×10.5	235	1pcs/bag
Drainage tubing	NPWT-03/200cm	Single use, Sterile	With blue pump and connectors	3 years	200 cm	120	1pcs/bag

Accessories of RH-1000 NPWT device



Canister holder



Canister bag 1000ml



Drainage tubes

- Instrument classification：Class I
- Scope of Application ：suitable for Acute & Chronic wounds, diabetic foot ulcers, venous ulcers, bedsores, neuropathic ulcers, donor sites, omphalocele, post-operative abdominoplasty ， etc.
- **It is combined with welsuc NPWT Device.**

Product name	Material code	Size	Futures	Air flow	Shelf life	Application	Weight(g)	Packaging Specifications
Canister holder	YLP01000Y	1000ml	Reusable, unsterile	<30L/min	3 years	Fixed the canister bag	230	20pcs/carton
Canister bag	YLP01000Y	1000ml	Disposable & unsterile with solidifier	<30L/min	3 years	Collet the excudate liquid	198	50pcs/carton
Drainage tube	YLP01000Y	180 cm	Disposable, sterile	<30L/min	3 years	Connected with canister bag and tube with pad	90	100pcs/carton

Negative Pressure Wound Therapy Dressing Kit

Black PU sponge with Silver avaialbel.



Each dressing kit including
1pcs black PU sponge
2pcs transparent film
1pcs dressing tubing set
1pcs I+Y connector
1pcs clamp

Black PU foam dressing



Each dressing kit including
1pcs white PVA sponge
2pcs transparent film
1pcs I+Y connector
1pcs clamp

PVA Foam Dressing

- Instrument classification：Class IIb
- Application：suitable for Acute & Chronic wounds, diabetic foot ulcers, venous ulcers, bedsores, neuropathic ulcers, donor sites, omphalocele, post-operative abdominoplasty，etc.

Product name	Material code	Dimensions (cm)	Material	Feature	Color	Aperture	Suitable for	Shelf life	sterile	Pack size
NPWT dressing kit	C.JE.WM.5053B	5×5×3	PU	Hydrophobicity	white/black	500~650μm	Chronic wound (diabetic foot, ulcer)	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.10073B	10×7×3	PU	Hydrophobicity	white/black	500~650μm	Chronic wound (diabetic foot, ulcer)	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.15103B	15×10×3	PU	Hydrophobicity	white/black	500~650μm	Chronic wound (diabetic foot, ulcer)	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.201003B	20×10×3	PU	Hydrophobicity	white/black	500~650μm	Chronic wound (diabetic foot, ulcer)	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.251503B	25×15×3	PU	Hydrophobicity	white/black	500~650μm	Chronic wound (diabetic foot, ulcer)	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.302003B	30×20×3	PU	Hydrophobicity	white/black	500~650μm	Chronic wound (diabetic foot, ulcer)	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.B10071A	10×7×1(or2or3)	PVA	Hydrophilicity	white	100~300μm	Burn, donor sites	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.B15101A	15×10×1(or2or3)	PVA	Hydrophilicity	white	100~300μm	Burn, donor sites	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.B20101A	20×10×1(or2or3)	PVA	Hydrophilicity	white	100~300μm	Burn, donor sites	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.B25151A	25×15×1(or2or3)	PVA	Hydrophilicity	white	100~300μm	Burn, donor sites	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.B30201A	30×20×1(or2or3)	PVA	Hydrophilicity	white	100~300μm	Burn, donor sites	2 years	Yes	1 set/bag

Remark: (1)This products are sterilized and disposable.
(2)The above specifications can be customized, if the quantity of each product is sufficient.
(3)Accessories of NPWT dressing kit:black/white foam, tubes With the pad, adhesive film, I/Y connector and clip.

Welsuc Negative Pressure Wound Therapy System

Model: RH-S01 & RH-S01I





- Instrument classification：Class II
- Scope of Application ：suitable for Acute & Chronic wounds, diabetic foot ulcers, venous ulcers, bedsores, neuropathic ulcers, donor sites, omphalocele, post-operative abdominoplasty，etc.

Advantages:

- ☐ Small & Portable.
- ☐ Three treatment modes
- ☐ Minor noise as working
- ☐ Extended working time with battery
- ☐ Smart operation
- ☐ Automatic instillation therapy
- ☐ Manage patients' therapy record

Specification of RH-S01 & RH-S01I dressing kit:

Product name	Material code	Dimensions (cm)	Model	Feature	Color	Time	Suitable for	Shelf life	sterile	Pack size
NPWT dressing kit	C.JE.WM.1020	10 ×20	A-1	Lock & solidifier liquid	white	2~7 days	low to moderately exuding wounds	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.1515	15×15	A-1	Lock & solidifier liquid	white	2~7 days	low to moderately exuding wounds	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.1520	15×20	A-1	Lock & solidifier liquid	white	2~7 days	low to moderately exuding wounds	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.1720	17×22	A-1	Lock & solidifier liquid	white	2~7 days	low to moderately exuding wounds	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.1030	10×30	A-1	Lock & solidifier liquid	white	2~7 days	low to moderately exuding wounds	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.1530	15×30	A-1	Lock & solidifier liquid	white	2~7 days	low to moderately exuding wounds	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.1725	17×30	A-1	Lock & solidifier liquid	white	2~7 days	low to moderately exuding wounds	2 years	Yes	1 set/bag
NPWT dressing kit	C.JE.WM.1040	10×40	A-1	Lock & solidifier liquid	white	2~7 days	low to moderately exuding wounds	2 years	Yes	1 set/bag

Remark: (1)This products are sterilized and disposable.
(2)The above specifications can be customized, if the quantity of each product is sufficient.
(3)Accessories of NPWT dressing kit:black/white foam, tubes With the pad, adhesive film, I/Y connector and clip.

RH 润虹医药 RAINHOME

Welsuc Negative Pressure Wound Therapy Dressing Kit

Clinical cases:

Application for the wounds as below:

- Acute or subacute wounds.
- Chronic wounds
- Dehisced wounds
- Pressure ulcers
- Diabetic/Neuropathic ulcers
- Venous insufficiency ulcers
- Traumatic wounds
- Partial thickness burns
- Flaps and grafts
- Closed surgical incisions













- Instrument classification：Class IIb for the device
- Negative pressure wound therapy (NPWT), also referred to as vacuum-assisted wound closure, is a therapeutic dressing system in which negative pressure is continuously or intermittently applied to the surface of a wound. The system includes dressings, a suction pump, tubing and a collection chamber.

润虹医药

Hydrocolloid Oil Gauze



Product name	Material code	Tape	Dimensions (cm)	Shelf life	Packaging Specifications
Hydrocolloid Oil Gauze	SJTC1010	Non-adhesive	10×10	3 years	1pcs/bag, 10bags/box
Hydrocolloid Oil Gauze	SJTC2010	Non-adhesive	20×10	3 years	1pcs/bag, 10bags/box

- Instrument classification：Class III
- Scope of Application ：Treatment of acute wounds (eg partial-thickness burns, donor site wounds, postoperative wounds and skin abrasions), chronic wounds (superficial ulcer and pressure ulcer care)
- Applicable departments ： Emergency department, ICU, general surgery, operating room, obstetrics and gynecology, burn department, endocrinology, neurosurgery, etc.

Advantages:

- Promote wounds healing, maintain moist environment for the wounds.
- Non-adhesive to the wound surface, pain relief and effectively protects the new tissues.
- Helps rebuild the dermis and epidermis.
- Safe and secure, no skin irritations, hypoallergenic,

Drug free, used in combination with other antibiotics.

Hydrocolloid Dressing



Product name	Material code	Tape	Dimensions (cm)	Shelf life	Packaging Specifications
Hydrocolloid Dressing	SJTD1010F	Adhesive	10×10	3 years	1pcs/bag, 10bags/box
Hydrocolloid Oil Gauze	SJTG1515F	Adhesive	15×15	3 years	1pcs/bag, 10bags/box

- Instrument classification：Class III
- Scope of Application ： Patients with low to moderate exudative wounds, including burns and donor site wounds; local ecchymosis, superficial ulcers; infusion phlebitis prevention and treatment and infusion needle fixation, etc.
- Applicable departments ： Emergency department, ICU, general surgery, operating room, obstetrics and gynecology, burn department, neurology, neurosurgery, etc.

Advantages:

- Absorb exudate, maintain moisture balance, and maintain the moist environment required for wound healing
- keep the wound surface moist and prevent the wound from scabbing
- Has a scab-dissolving effect on scabbed wounds and reduces scarring
- It can prevent the product from sticking to the wound and reduce secondary damage.
- Occlusive dressing with antibacterial effects
- Breathable film, less irritating to skin, low sensitization, good adhesion.

Alginate Dressing



Product name	Material code	Tape	Dimensions (cm)	Shelf life	Packaging Specifications
Alginate Dressing	ZSF0706	Non-adhesive	7×6	2 years	1pcs/bag, 10bags/box
Alginate Dressing	ZSF1010	Non-adhesive	10×10	2 years	1pcs/bag, 10bags/box
Alginate Dressing	ZSN1010	Adhesive	10×10	2 years	1pcs/bag, 10bags/box
Alginate Dressing	ZSN1510	Adhesive	15×10	2 years	1pcs/bag, 10bags/box
Alginate Dressing	ZSF1010E	Non-adhesive	10×10	2 years	1pcs/bag, 10bags/box
Alginate Dressing	ZSN1010L	Adhesive	10×10	2 years	1pcs/bag, 10bags/box
Alginate Dressing	ZSN1510M	Adhesive	15×10	2 years	1pcs/bag, 10bags/box

- Instrument classification：Class III
- Scope of Application ：Suitable for covering on the wounds, absorbing the wound exudate. - Suitable for wounds with cavities. After filling, it absorbs exudate and provides a moist environment for wounds healing.
- Applicable departments ： General surgery, anorectal department, obstetrics and gynecology, burn and plastic surgery, ENT, endocrinology, neurosurgery, nursing ward, etc.

Advantages:

- **100% pure alginate without others chemical ingredients added and Non-toxic.**
- Does not cause skin allergies, does not cause infection
- Non-sticking to the skin, easy to remove, less bleeding rate during dressing change, relief pain.
- Soft material, form a soft gel & comfortable use.
- Promotes wounds healing, no stimulation & less scar formation.

RAINHOME 润虹医药

Gelling Fiber Dressing



Product name	Material code	Tape	Dimensions（cm）	Shelf life	Packaging Specifications
Gelling Fiber Dressing	QSA0505	Non-adhesive	5×5	2 years	1pcs/bag, 10bags/box
Gelling Fiber Dressing	QSC1010	Non-adhesive	10×10	2 years	1pcs/bag, 10bags/box

- Instrument classification：Class III
- Scope of Application： It is used for acute and chronic wounds with moderate to severe exudation, and can also be used for dry or mildly exuding wounds after pre-wetting treatment.
- Applicable departments ： Neurosurgery, Cardiac Surgery, Thoracic Surgery, General Surgery, Obstetrics and Gynecology, Urology, Orthopedics, ENT, etc.

Advantages:

- Composed of natural fibers (cellulose derivatives), with strong exudate management ability
- Soft application in gel state, reducing dead space and reducing bacteria living space
- Reduces pain and does not damage tendons, wound granulation tissue, or healthy tissue surrounding the wound.
- After absorbing the exudate, it becomes transparent and gelatinized, with a high degree of transparency, and the wound healing can be observed.

RAINHOME 润虹医药

Foam Wound Dressing (with silver available)



Product name	Material code	Tape	Dimensions (cm)	Shelf life	Packaging Specifications
Foam Dressing	PMB1010	Non-adhesive	10×10	2 years	1pcs/bag, 10bags/box
Foam Dressing	PMB1515	Non-adhesive	15×15	2 years	1pcs/bag, 10bags/box
Foam Dressing	PMC1010N	Adhesive	10×10	2 years	1pcs/bag, 10bags/box
Foam Dressing	PMC1515N	Adhesive	15×15	2 years	1pcs/bag, 10bags/box

- Instrument classification：Class III

- Scope of Application ： All kinds of moderate to severe exuding wounds, dispersing pressure, preventing pressure ulcers, combined with alginate in donor sites, especially for people with fragile skin.

- Applicable departments ：ICU/NICU, Orthopedics, burns, anorectal, endocrinology, hand surgery, general surgery, wound clinic, interventional department, pediatrics, geriatrics, etc.

Advantages:

- Absorbs exudate, promotes tissue growth, thickens the foam layer deeply locks the liquid to keep the skin around the wound dry and not soaked.
- Antibacterial & breathable, prevents wounds infection.
- Relief pain as changing, Non-irritating and promotes wounds healing. Reduces the secondary injury and mechanical damage of the wound.
- Constant temperature, humidity, hypoxia and slightly acid environment that can promote release of growth factors and the movement of wound cells fast.

RAINHOME 润虹医药

Chitosan Wound Dressing



- Instrument classification：Class III
- Scope of Application ： used for pressure ulcers, Graft & donor site, arterial & diabetic ulcers, leg ulcers, postoperative surgical wounds, burns, cavity wounds.
- Applicable departments ：Orthopedics, burns, anorectal, endocrinology, hand surgery, general surgery, wound clinic, interventional department, pediatrics, geriatrics, etc.

Product name	Material code	Tape	Dimensions (cm)	Shelf life	Packaging Specifications
Chitosan Wound Dressing	KJF0610	Non-adhesive	6×10	2 years	1pcs/bag, 10bags/box

Advantages:

- ☐ Natural material & promote wounds healing.
- ☐ Moisturizing and breathable.
- ☐ Antibacterial & reduce scar formation.
- ☐ Stops bleeding and pain.
- ☐ Prevent wounds adhesion
- ☐ Great safety & green product.

RAINHOME 润虹医药

Medical Scar Paster



- Instrument classification：Class II
- Scope of Application ： Suitable for preventing and inhibiting scarring.
- Applicable departments ： Plastic surgery, burns, dermatology, obstetrics and gynecology, cardiothoracic surgery, orthopedics, breast surgery, pediatric surgery, ENT, general surgery, etc.

Product name	Material code	Tape	Dimensions (cm)	Shelf life	Packaging Specifications
Medical Scar Paster	BH00123L	Non-adhesive	12×3	3 years	1pcs/bag, 10bags/box
Medical Scar Paster	BH00128	Non-adhesive	12×8	3 years	1pcs/bag, 10bags/box

Advantages:

- ☐ Great safety, non-toxic and non-irritating.
- ☐ Good Breathable.
- ☐ Soft texture, easy to cut.
- ☐ Good skin affinity and hydrophilicity.
- ☐ Long lasting effect.

Medical Functional Dressing

Medical Functional Dressing



Specifications: 15ml/bottle, 20ml/bottle, 30ml/bottle, 50ml/bottle/box.

Main components: Chitosan

- Instrument classification：Class II
- Scope of Application ：Suitable for prevention of the wounds infection of skin and mucous membrane, and promotes the wounds healing.
- Applicable departments ： Emergency clinic, general surgery, obstetrics and gynecology, ENT, stomatology, dermatology, burn and plastic surgery, anorectal, oncology, etc.

Features:

- Antibacterial, anti-inflammatory, prevent infection.
- Stop bleed, pain relief, Moisturizing, Breathable.
- Repair scars：Promotes tissues repair and regeneration, inhibit and repair scars.

Saline Nasal Spray





- Instrument classification：Class II
- Scope of Application ： Suitable for cleaning the nasal cavity and keeping the nasal cavity moist.
- Applicable departments ： Otolaryngology, pediatrics, neonatology, respiratory medicine, infectious diseases, oncology (after nasopharyngeal radiotherapy), etc..

Product name	Material code	Tape	Dimensions (ml)	Shelf life	Packaging Specifications
Saline Nasal Spray	YS00030H	E	30	3 years	180pcs/carton
Saline Nasal Spray	YS00060H	M	60	3 years	180pcs/carton
Saline Nasal Spray	YS00050H	K	50	3 years	100pcs/carton
Saline Nasal Spray	YS00030	/	30	3 years	300pcs/carton
Saline Nasal Spray	YS00060	/	60	3 years	180pcs/carton

Advantages:

- Cleanse and keep the nasal cavity moist,.
- Quick relief due to allergic rhinitis and flu.
- More economical, safe and portable.

Marketed 18 years , Sold 1500W+pcs

润虹医药

Medical ultrasound gel

Medical ultrasound gel for cavity and Medical ultrasound gel for disinfection



- Instrument classification：Class II
- Scope of Application ： Suitable for ultrasound examination before and after surgery in a sterile environment;Ultrasound examination of cavity mucosa.
- Applicable departments ： Obstetrics and gynecology, gastroenterology, urology, neurology, neonatology, etc.

Product name	Material code	Tape	Dimensions (g)	Shelf life	Packaging Specifications
Medical ultrasound gel for cavity	OHWJ0012	Sterilized	12	2 years	400pcs/carton
Medical ultrasound gel for cavity	OHWJ0020	Sterilized	20	2 years	400pcs/carton
Medical ultrasound gel	OHXD0012	disinfection	12	2 years	400pcs/carton
Medical ultrasound gel	OHXD0020	disinfection	20	2 years	400pcs/carton
Medical ultrasound gel	OH00250	disinfection	250	2 years	60pcs/carton

Advantages:

- ☐ Good sound transmission function and coupling function make the ultrasound image clearer.
- ☐ Strong sterilization function.
- ☐ More Safe and non-irritating.
- ☐ No corrosion to the probe

Marketed 18 years , Sold 1500W+pcs

润虹医药

Thank You

Guangzhou Rainhome Pharm & Tech Co., Ltd
Contact person: Jason Huang
Email: ds03@rhkj.com.cn
Tel: +86-13632492562 (Wechat/WhatsApp)